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User Manual

Thank you for choosing our product. We hope you enjoy using the appliance.

Contents

General safety Instructions	35
Symbols in these Instructions for Use	35
Children.....	35
Special Safety Instructions For This Machine	35
Handling Batteries	35
Inserting the Batteries	36
Location of Controls	36
Technical Specifications	36
Information on Measuring Body Fat	37
Body Fat.....	37
The Importance of Water for our Health.....	37
Guideline Figures.....	37
Use	38
Weighing Function	38
Saving Personal Data	38
Measurement and Analysis.....	38
Display Error Messages	39
Cleaning	39
Disposal - Meaning of the "Dustbin" Symbol	39
Disposal of batteries.....	39

General safety Instructions

Read the operating instructions carefully before putting the appliance into operation and keep the instructions including the warranty, the receipt and, if possible, the box with the internal packing. If you give this device to other people, please also pass on the operating instructions.

- The device is designed exclusively for private and not for commercial use.
- Only use the device for the intended purpose.
- Do not use a damaged device.
- Keep away from heat, direct sunlight, damp (never immerse in liquids) and sharp edges.
- Do not try to repair the appliance on your own. Always contact an authorized technician.
- Pay careful attention to the following "Special Safety Instructions".

Symbols in these Instructions for Use

Important information for your safety is specially marked. It is essential to comply with these instructions in order to avoid accidents and prevent damage to the machine:

WARNING:

This warns you of dangers to your health and indicates possible injury risks.

CAUTION:

This refers to possible hazards to the machine or other objects.

NOTE:

This highlights tips and information.

Children

- In order to ensure your children's safety, please keep all packaging (plastic bags, boxes, polystyrene etc.) out of their reach.

WARNING!

Caution! Do not allow small children to play with the foil as there is a **danger of suffocation!**

- Children should be supervised at all times in order to ensure that they do not play with the device.

Special Safety Instructions For This Machine

WARNING:

If the glass surface is moist or covered with water droplets, please wipe the surface dry as there is a **danger of slipping**.

CAUTION:

- Handle the balance with care, as it is a precise measuring instrument. Do not drop it or jump on it.
- Do not load the appliance with more than 150 kg (330.7 lb, 23.62 st) to avoid damaging the device.
- If the scales are carried or not going to be used, please store them only in a horizontal position.
- The body fat analysis scale is not suitable for the following persons:
 - Children under 10 years and adults over 85 years
 - Persons with symptoms of fever, oedema, or osteoporosis
 - Persons undergoing treatment for dialysis
 - Persons who use a pacemaker
 - Persons who take cardiovascular medication
 - Pregnant women
 - Sportspersons who perform over 5 hours of intensive sporting activities each day
 - Persons with a heartbeat lower than 60 beats/minute

Handling batteries

WARNING:

- Do not expose the batteries to high temperatures or direct sunlight. Never throw batteries into the fire. There is a **danger of explosion!**
- Keep batteries away from children. They are **not** a toy!
- Do not open the batteries by force.

36 English

⚠ WARNING:

- Avoid contact with metallic objects. (Rings, nails, screws et cetera.) There is a **danger of short-circuiting!**
- As a result of a short-circuit batteries may heat up considerably or even catch fire. This may result in burns.
- For your own safety the battery poles should be covered with adhesive strips during transport.
- If a battery leaks do not rub the liquid into your eyes or mucous membranes. If hands come into contact with the liquid they should be washed thoroughly and eyes rinsed out with clear water; if symptoms persist please consult a doctor.

Inserting the batteries

📌 NOTE:

Where required, the supplied batteries have been protected for transport with plastic foil.

1. Remove this foil before inserting the batteries into the appliance.
2. Open the battery compartment on the back.
3. Insert 2 batteries of type AAA/R03. Please ensure that the polarity is correct (see the loor of the battery compartment)!
4. Close the battery compartment again.
5. If the device is not used for a prolonged period of time, please take the batteries out in order to prevent leakage of battery acid.
6. If the message "Lo" appears on the display, replace the batteries.

⚠ CAUTION:

- Different battery types or new and used batteries must not be used together.
- Batteries are not to be disposed of together with domestic wastes. Dispose of the used batteries at the relevant collection bin or at your dealer's.

Location of controls

- 1 Contact Surfaces
- 2 Display
- 3 DOWN-button ▼
- 4 UP-button ▲
- 5 SET-button

Display symbol	Meaning
kg (lb/st)	Body weight
USER	User memory space
	Male

Display symbol	Meaning
	Female
cm (Ft:in)	Body height
Age	Age
	Unathletic type
	Average type
	Athletic Type
	Body fat display
	Body water display
	Muscle mass display
	Bone weight display
%	Percentage – body fat / body water / muscle mass / bone weight
	Analysis bar: Extremely underweight – underweight – normal weight – overweight – obese
KCAL	Daily kilocalorie requirement
Increase	Weight gain in kg (lb/st) compared with previous measurement
Decrease	Weight loss in kg (lb/st) compared with previous measurement

Technical Specifications

Model:PW 5571 FA
 Power supply:..... 2 Batteries 1.5 V Type "AAA" / "R03"
 Net weight of device 2.5 kg
 Memory capacity:.....10
 Mass alignment:from 75 to 225 cm (2'5" – 7'4½" Ft:in)
 Age:..... from 10 to 85 years
 Minimum weight:from 2.5 kg (5.5 lb/0: 4 st)
 Maximum weight: 150 kg (330.7 lb/23: 6 st)
 Weight measurement unit: 100 g increments
 Body fat measurement unit:..... 0.1 %
 Body fat measuring range:.....4.0 % - 50.0 %
 Body water measuring range:35.0 % - 85.0 %
 Muscle mass measuring range:..... 10.0 % - 50.0 %
 Subject to technical changes without prior notice!

This device has been tested according to all relevant current CE guidelines, such as electromagnetic compatibility and low voltage directives, and has been constructed in accordance with the latest safety regulations.

Information on Measuring Body Fat

WARNING:

Do not use the scales if you use a pacemaker or another medical implant with electronic components. The function of the implant may be impaired by the measurement.

NOTE:

The device can only function correctly if you stand on the contact surfaces with bare, dry feet!

The main principle of this device is based on the measurement of the electric impedance within the human body. In order to determine the actual body data, the device leads a weak, unnoticeable electric signal through your body.

The signal is generally safe! However, please read our special safety instructions!

The method of measurement is referred to as "bioelectric impedance analysis" (BIA), a factor which is based on the relationship between body weight and associated body water percentages as well as on other biological data (age, gender, body weight).

There are heavy persons with a well-balanced percentage between fat and muscle mass, and there are seemingly slim persons with significantly high body fat values. Fat is a factor that puts serious strain on the body, in particular the cardiovascular system. Therefore, controlling both values (fat and weight) is a prerequisite for health and fitness.

Body Fat

Body fat is an important component of the body. From protecting the joints as a vitamin repository to regulating the body temperature, fat plays a vital role in the body. Thus the goal is not to radically reduce the proportion of body fat, but rather to achieve and maintain a healthy balance between muscle and fat.

The Importance of Water for our Health

Approximately fifty-five to sixty percent of the human body consist of water, depending on age and gender in varying amounts.

Water has to fulfill several tasks within the body:

- It is the building block of our cells. That means that all the cells in our body, whether skin cells, gland cells, muscle cells, brain cells or others, are functional only if they contain enough water.
- It is a solvent for important substances in our body.
- It is a means of conveyance for nutrients, the body's own substances and metabolic products.

Reference values for your fitness profile can be seen in the benchmark table. It shows the dependency of body fat data on age and gender. Always pay attention to your body's water balance, which, as well as the fat level, undergoes strong fluctuations throughout the day. Therefore, to obtain consistent data, always measure at the same time of day.

The calculated results of the body fat and water level should only be considered as reference values. They are not recommended to be used for medical purposes. Please always consult a doctor if you have questions about a medical treatment or a diet in order to achieve a normal body weight.

Guideline values for body fat and body water content as a percentage (approximate figures)

Age (Years)	Women		Men		
	Body fat	Body water	Body fat	Body water	
< = 30	4.0 – 20.5	66.0 – 54.7	4.0 – 15.5	66.0 – 58.1	Underweight
	20.6 – 25.0	54.6 – 51.6	15.6 – 20.0	58.0 – 55.0	Normal
	25.1 – 30.5	51.5 – 47.8	20.1 – 24.5	54.9 – 51.9	Overweight
	30.6 – 50.0	47.7 – 34.4	24.6 – 50.0	51.8 – 34.4	Obese
> 30	4.0 – 25.0	66.0 – 51.6	4.0 – 19.5	66.0 – 55.4	Underweight
	25.1 – 30.0	51.5 – 48.1	19.6 – 24.0	55.3 – 52.3	Normal
	30.1 – 35.0	48.0 – 44.7	24.1 – 28.5	52.2 – 49.2	Overweight
	35.1 – 50.0	44.6 – 34.4	28.6 – 50.0	49.1 – 34.4	Obese

38 English

Use

The device has a very sensitive sensor. If you want to adjust the settings, simply touch the buttons gently with your finger. If the device does not respond to touch, clean the glass surface with a soft, dry cloth.

i NOTE:

The scales are factory set with the units "kg" and "cm".

To change the settings to "lb" or "st" as well as "Ft.in", proceed as follows:

- Switch the scales on by applying gentle pressure on the glass surface with your foot or toe.
- Press the UP button ▲, until the required unit of weight is displayed. The unit for the height also automatically changes.

Weighing Function

- Place the scale on a hard, flat surface (avoid carpets). An uneven floor contributes to inaccurate measurements.
- Stand on the scales with both feet. The device switches on automatically.
- Distribute your weight equally and wait for the scale to determine your weight.
- The digits on the display will flash before your exact weight will be shown.
- In order to extend the battery life, the device will automatically turn off 10 seconds after you have stepped off the scale.

Saving Personal Data

i NOTE:

- The scale has default values that you can change, as described below.
- Carry out the following steps in sequence and within a few seconds. Otherwise, the device will automatically cancel input mode.
- Press the button SET. The memory space will flash on the display (USER). Confirm the setting by pressing SET or select another memory space by pressing the UP / DOWN button. Then press SET to go to the next setting.
- The symbol for "male" or "female" will flash. Confirm the setting by pressing SET or use the UP / DOWN button to choose the symbol for your gender. Then press SET to go to the next setting.
- The height setting will flash. Select your height by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.
- The setting AGE will flash. Select your age by pressing UP or DOWN button. Holding down one of the buttons will accelerate the process. Confirm your entry by pressing SET.

- The default symbol for the "unathletic" type flashes. Which type do you consider yourself to be? Select your type using the UP or DOWN button: "unathletic", "average" or "athletic". Confirm your selection with SET.
- If only the stored position and the display 0.0 kg (0.0 lb/0.0 st) is still shown, input mode for this stored position is completed.

You now have the following options:

- You can change the units for weight and type. Press the UP button ▲, until the desired unit of weight is displayed. The unit for the height automatically changes.
- Press the SET button to see the selected settings again.
- Allocate further memory spaces. There are 10 spaces available in total. Press SET and change the space by pressing the UP or DOWN button. Proceed further as described above.
- If you want to make changes, press SET (USER flashes). Select the stored position using the UP or DOWN button. Proceed further as described above.
- You can immediately start measuring.

Measurement and Analysis

Record your personal data in the memory.

i NOTE:

- The device can only function correctly if you stand on the device with bare, dry feet!
- If necessary wipe the contact surfaces to improve results.
- Due to the significant fluctuation of body fat, it is difficult to attain a standard value. On average the optimal testing time is between 6 p.m. and 8 p.m, since the body fat remains very consistent during this period of time.
- Press the SET button. The preset stored position data appears on the display (USER flashes). Use the UP / DOWN button to select the stored position at which your data is saved.
- Wait for the USER display and 0.0 kg (0.0 lb/st) to stop flashing.
- Stand on the contact surfaces (1) of the scales with both feet.
- Firstly your weight will flash on the display.
- After a few seconds, the following results will be shown:
 1. - the body fat content in % (☞ symbol)
- the body water content in % (☞ symbol)
- the results analysis in the form of a bar
 2. - the muscle mass content in % (☞ symbol)
- the bone mass content in % (☞ symbol)
- the calorie requirement per day
 3. - The preset gender
- The preset type
- The preset height
- The preset age (AGE)

- The weight gain (INCREASE) / the weight loss (DNCREASE) in kg (lb/st) in comparison to the previous measurement
- Afterwards the determined data will be shown repeatedly. Then the scale will automatically turn off.

Display error Messages

i NOTE:

If an error message appears on the display, wait until the device is turned off. Turn it on again to restart the operation. If the error display has not disappeared, take the batteries out of the battery compartment for 10 seconds.

Display symbol	Meaning
8888	Self test (when switching on) to check that the numeric display functions properly
----	The scales stand faultily
Lo	Battery is low
Err	The percentage of body fat content lies above or below the limits of the measuring scale
EEEE	The weight lies above the allowed maximum weight.

Cleaning

⚠ CAUTION:

- Do not immerse the balance in water: the electronic components can be damaged.
- Do not use a wire brush or any abrasive items.
- Do not use any acidic or abrasive detergents.
- Clean the device with a slightly moist cloth. Make sure that no humidity or waters seeps into the device!
- Use a commercially sold detergent or soap for cleaning.
- Dry the device well after cleaning with a dry and soft cloth.

Disposal - Meaning of the "Dustbin" Symbol

Protect our environment: do not dispose of electrical equipment in the domestic waste.

Please return any electrical equipment that you will no longer use to the collection points provided for their disposal.

This helps avoid the potential effects of incorrect disposal on the environment and human health.

This will contribute to the recycling and other forms of reutilization of electrical and electronic equipment.

Information concerning where the equipment can be disposed of can be obtained from your local authority.

Disposal of batteries

Your appliance contains batteries. If these batteries are spent, follow the instructions below:

⚠ CAUTION:

Batteries/rechargeable batteries must not be disposed of with household waste!

- The batteries must be removed from the appliance.
- Take spent batteries to the appropriate collection point or to a dealer.
- Your town or local authority can provide information about public collection points.

This symbol can be found on batteries/rechargeable batteries which contain hazardous substances:

- Pb = contains lead
- Cd = contains cadmium
- Hg = contains mercury
- Li = contains lithium

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